






Oxivir[®] Tb Wipes

Hospital Grade One-Step Cleaner Disinfectant

	<ol style="list-style-type: none"> 1. Pull towelette from the canister dispenser 2. Wipe the surface 	<p>TGA REGISTERED – AUST R 164850</p> <p>Proven to kill:</p> <p>Viruses:</p> <ul style="list-style-type: none"> • H1N1 Influenza A (Swine Flu) • Hepatitis B • Hepatitis C* • Herpes simplex virus type1 & type 2 • HIV-1 (AIDS virus) • Human Coronavirus • Influenza A • Norovirus • Polio virus • Rhinovirus • Rotavirus • Avian Influenza A • Feline calicivirus <p>Bacteria:</p> <ul style="list-style-type: none"> • Acinetobacter baumannii • Escherichia coli⁺ • Escherichia coli 0157:H7 • Klebsiella pneumoniae • Proteus vulgaris⁺ • Pseudomonas aeruginosa • Salmonella enterica • Shigella Dysenteriae • Staphylococcus aureus <p>Antibiotic- Resistant Bacteria:</p> <ul style="list-style-type: none"> • Enterococcus faecalis (VRE) • Staphylococcus aureus (MRSA) <p>Fungi:</p> <ul style="list-style-type: none"> • Trichophyton mentagrophytes⁺ <p>Tb:</p> <ul style="list-style-type: none"> • Mycobacterium bovis* <p>1 minute contact time for viruses & bacteria * 5 minutes contact time ⁺ 10 minutes contact time</p>
	<ol style="list-style-type: none"> 1. All surfaces must remain visibly wet for 1 minute to have a virucidal and bactericidal effect 2. Allow to air dry 3. Dispose wipe after use 4. Do not flush in a toilet 	
	<ol style="list-style-type: none"> 1. Ideal for keyboards, telephones, patient areas and high-touch surfaces 2. If disinfecting food contact surfaces, rinse with water. 	
<p><u>SUPPLIES/TOOLS:</u></p> <ul style="list-style-type: none"> • Gloves should be considered, as prolonged hand contact will dry skin 	<p><u>SAFETY:</u></p> <ul style="list-style-type: none"> • Oxivir Tb has a low concentration of active (AHP) 	<p><u>CONTROL POINTS:</u></p> <ul style="list-style-type: none"> • Do not mix with other detergents or other chemicals • Not for use on glasses, dishes and utensils or therapeutic devices. • Not recommended for use on brass, copper or marble • Safe on carpets and fabrics • Safe on stainless steel • Biodegrades to OXYGEN and WATER