






# Oxivir<sup>®</sup> Tb Wipes

## Hospital Grade One-Step Cleaner Disinfectant

	<ol style="list-style-type: none"> <li>1. Pre-clean heavily soiled areas</li> <li>2. Pull towelette from the dispenser (canister)</li> <li>3. Wipe the surface</li> </ol>	<p><b>TGA EPA Registration No. 70627-60</b></p> <p>Tested &amp; proven effective against a wide spectrum of pathogens:</p> <p><b>Viruses:</b></p> <ul style="list-style-type: none"> <li>• Hepatitis B</li> <li>• Hepatitis C</li> <li>• Herpes Simplex Virus, Type 1 &amp; 2</li> <li>• HIV-1 (AIDS virus)</li> <li>• Human Coronavirus</li> <li>• Influenza A</li> <li>• Norovirus</li> <li>• Poliovirus</li> <li>• Rotavirus</li> <li>• Adenovirus</li> <li>• Rhinovirus Type 14</li> <li>• H1N1 Influenza A (Swine flu)</li> <li>• Avian influenza A</li> <li>• Feline Calicivirus</li> </ul> <p><b>Bacteria:</b></p> <ul style="list-style-type: none"> <li>• Pseudomonas aeruginosa</li> <li>• Salmonella enterica</li> <li>• Staphylococcus aureus</li> <li>• Acinetobacter baumannii</li> <li>• Escherichia coli O157:H7</li> <li>• Klebsiella pneumoniae</li> <li>• Escherichia coli with extended beta-lactamase resistance (ESBL)</li> <li>• Shigella Dysenteriae</li> </ul> <p><b>Antibiotic- Resistant Bacteria:</b></p> <ul style="list-style-type: none"> <li>• Enterococcus faecalis (VRE)</li> <li>• Staphylococcus aureus (MRSA)</li> <li>• Staphylococcus aureus (CA-MRSA)</li> <li>• Staphylococcus aureus (CA-MRSA)</li> </ul> <p><b>Fungi:</b> Trichophyton mentagrophytes</p> <p><b>Tb:</b> Mycobacterium bovis</p>
	<ol style="list-style-type: none"> <li>1. All surfaces must remain visibly wet for 1 minute. for viruses and bacteria. Allow 5 minutes for Tb and 10 minutes for fungi.</li> <li>2. Allow to air dry</li> <li>3. Dispose wipe after use.</li> <li>4. Do not flush in a toilet</li> </ol>	
	<ol style="list-style-type: none"> <li>1. Ideal for keyboards, telephones, patient areas and high-touch surfaces</li> <li>2. If disinfecting food contact surfaces, rinse with water.</li> </ol>	
<p><b><u>SUPPLIES/TOOLS:</u></b></p> <ul style="list-style-type: none"> <li>• Gloves should be considered, as prolonged hand contact will dry skin</li> </ul>	<p><b><u>SAFETY:</u></b></p> <p><b>No PPE Gear Needed.</b></p> <ul style="list-style-type: none"> <li>• Oxivir Tb has a low concentration of active</li> <li>• Wear utility gloves and other PPE appropriate to the task such as site decontamination after spills of blood or other potentially infectious materials (NHMRC Guidelines).</li> </ul>	<p><b><u>CONTROL POINTS:</u></b></p> <ul style="list-style-type: none"> <li>• Do not mix with other detergents or other chemicals</li> <li>• Not for use on glasses, dishes and utensils or therapeutic devices.</li> <li>• Not recommended for use on brass, copper or marble</li> <li>• Safe on carpets, fabrics &amp; clothes</li> <li>• Safe on stainless steel</li> <li>• Post use the active ingredient breaks down to only oxygen and water</li> </ul>