






Oxivir[®] Tb

Hospital Grade One-Step Cleaner Disinfectant

| | | |
|--|---|--|
|  | <ol style="list-style-type: none"> 1. Apply Oxivir Tb undiluted to surfaces by spray, cloth or disposable wipe. | <p>TGA REGISTERED – AUST R 165058</p> <p>Tested & proven effective against a wide spectrum of pathogens:</p> <p>Viruses:</p> <ul style="list-style-type: none"> • H1N1 Influenza A (Swine Flu) • Hepatitis B • Hepatitis C* • Herpes simplex virus type1 & type 2 • HIV-1 (AIDS virus) • Human Coronavirus • Influenza A • Norovirus • Rhinovirus • Rotavirus • Polio virus • Avian Influenza A • Feline calicivirus <p>Bacteria:</p> <ul style="list-style-type: none"> • Acinetobacter baumannii • Escherichia coli⁺ • Escherichia coli 0157:H7 • Klebsiella pneumoniae • Proteus vulgaris⁺ • Pseudomonas aeruginosa • Salmonella enterica⁺ • Shigella Dysenteriae • Staphylococcus aureus <p>Antibiotic- Resistant Bacteria:</p> <ul style="list-style-type: none"> • Enterococcus faecalis (VRE) • Staphylococcus aureus (MRSA) <p>Fungi:</p> <ul style="list-style-type: none"> • Trichophyton mentagrophytes⁺ <p>Tb:</p> <ul style="list-style-type: none"> • Mycobacterium bovis* <p><small>1 minute contact time for viruses & bacteria * 5 minutes contact time * 10 minutes contact time</small></p> |
|  | <ol style="list-style-type: none"> 1. All surfaces must remain visibly wet for 1 minute for viruses and bacteria (refer to the PSS or label for specific organisms and contact time) 2. Allow to air dry. 3. If concerned about streaking, wipe with a microfibre cloth. | |
|  | <ol style="list-style-type: none"> 1. If disinfecting food contact surfaces, rinse with water. 2. If Clostridium Difficile outbreak is identified, use microfibre . Or use Sporicide Plus. | |
| <p>SUPPLIES/TOOLS:</p> <ul style="list-style-type: none"> • Ultra-microfibre cloths • Disposable dry towel • Gloves should be considered, as prolonged hand contact will dry skin. | <p>SAFETY:</p> <p>No PPE Gear Needed.</p> <ul style="list-style-type: none"> • Oxivir Tb has a low concentration of active • Wear utility gloves and other PPE appropriate to the task such as site decontamination after spills of blood or other potentially infectious materials (NHMRC Guidelines). | <p>CONTROL POINTS:</p> <ul style="list-style-type: none"> • Do not mix with other detergents or other chemicals • Not for use on glasses, dishes and utensils or therapeutic devices. • Not recommended for use on brass, copper or marble • Safe on carpets, fabrics & clothes • Safe on stainless steel • Post use the active ingredient breaks down to only oxygen and water |